



www.diadent.co.kr



2014 / 2015
DiaDent Group International

Passion for Innovation ●●●

DiaDent's Network

DiaDent Group International(Korea & International)

No. 626, Yeonje-Ri, Gangoe-Myeon, Cheongwon-Gun, Chungcheong Buk Do, 363-951, South Korea
Tel)82-43-266-2315 Fax)82-43-262-8658
E-mail: diadent@diadent.co.kr

DiaDent Group International Inc. (Canada & USA)

11-3871 North Fraser Way, Burnaby B.C, Canada V5J 5G6
Tel)1-604-451-8851 Fax)1-604-451-8865
E-mail: diadent@diadent.com

DiaDent Europe B.V (Europe)

Antennestraat 70, 1322 AS Almere, the Netherlands
Tel)31-36-549-8607 Fax)31-36-536-7317
E-mail: diadent@diadent europe.com

REP

2013/11/20 [REV. 1]

©Copyright 2014 DiaDent Group International

www.diadent.co.kr

Welcome to DiaDent

Passion for Innovation!

DiaDent is one of the leading Endodontic supply companies in the dental industry since 1985. For more than 29 years, DiaDent has been dedicated to making a better and high-quality dental products through our passion for the Innovation and leads the global market in high-advanced Endodontic items manufacturing and supplying.

From its inception as exporting Gutta Percha Points and Paper Points, DiaDent has grown to become a global leader, specializing in Endodontic devices and materials. DiaDent's innovative and reliable items, services and processes are recognized as the world's brand, today DiaDent has been exporting them to over 900 dealers in 117 countries all around the globe. DiaDent is taking the world in different systems with the highest standards so that DiaDent can more effectively contribute to create new markets and continue to make a responsible approach to the global business. All of DiaDent's products are inspected by a strict Quality Management System following the ISO 9001 & ISO 13485 standards.

Always striving to be better and innovative, DiaDent puts in our utmost efforts to strengthen R&D part to develop brand-new dental devices and items. Our products from DiaDent's R&D part lead us to the global market in high-tech dental manufacturing. Through innovative, reliable products and services and a responsible approach to the business and collaboration with our partners and customers, DiaDent will continue to see every challenge as an opportunity and we will do the best to be perfectly positioned as one of the world's recognized leaders in the dental industry.

Everyone here at DiaDent would like to thank you, our valued customers and business partners, for your continued support and loyalty. We wish everyone another happy and prosperous year ahead.

Sincerely,

DiaDent Group International

DiaDent's Network

DiaDent Group International (Korea & International)

No. 626, Yeonje-Ri, Gangoe-Myeon, Cheongwon-Gun, Chungcheong Buk Do, 363-951, South Korea
Tel) 82-43-266-2315 Fax)82-43-262-8658
E-mail: diadent@diadent.co.kr

DiaDent Group International Inc. (Canada & USA)

11-3871 North Fraser Way, Burnaby B.C, Canada V5J 5G6
Tel)1-604-451-8851 Fax)1-604-451-8865
E-mail: diadent@diadent.com

DiaDent Europe B.V (Europe)

Antennestraat 70, 1322 AS Almere, the Netherlands
Tel)31-36-549-8607 Fax)31-36-536-7317
E-mail: diadent@diadent europe.com





DiaDent's innovative and reliable items, services and processes are recognized as the world's brand, today DiaDent has been exporting them to over 900 dealers in 117 countries all around the globe.

ALBANIA	MALAYSIA
ARMENIA	MALTA
ARGENTINA	MAURITIUS
AUSTRALIA	MEXICO
AUSTRIA	MOLDOVA
AZERBAIJAN	MONGOLIA
BAHRAIN	MONTENEGRO
BANGLADESH	MOROCCO
BELARUS	MYANMAR
BELGIUM	NETHERLANDS
BOLIVIA	NEW ZEALAND
BOSNIA & HERZEGOVINA	NEPAL
BRAZIL	NIGERIA
BRUNEI	NORWAY
BULGARIA	OMAN
CANADA	PAKISTAN
CAMBODIA	PALESTINE
COSTA RICA	PANAMA
CHILE	PAPUA NEW GUINEA
CHINA	PARAGUAY
COLOMBIA	PERU
CROATIA	PHILIPPINES
CYPRUS	POLAND
CZECH REPUBLIC	POLYNESIA
DENMARK	PORTUGAL
DOMINICAN REPUBLIC	QATAR
ECUADOR	REUNION ISLAND
EGYPT	ROMANIA
EL SALVADOR	RUSSIA
ESTONIA	SAUDI ARABIA
FINLAND	SENEGAL
FRANCE	SERBIA
GEORGIA	SINGAPORE
GERMANY	SLOVENIA
GHANA	SOUTH AFRICA
GREECE	SPAIN
GUATEMALA	SRI LANKA
HONG KONG	SUDAN
HUNGARY	SWEDEN
INDIA	SWITZERLAND
INDONESIA	SYRIA
IRAN	TAIWAN
IRAQ	THAILAND
IRELAND	TRINIDAD & TOBAGO
ISRAEL	TUNISIA
ITALY	TURKEY
JAMAICA	UAE
JAPAN	UGANDA
JORDAN	UNITED KINGDOM
KAZAKHSTAN	UKRAINE
KENYA	U.S.A
KOREA	URUGUAY
KOSOVO	UZBEKISTAN
KUWAIT	VENEZUELA
LATVIA	VIETNAM
LEBANON	YEMEN
LIBYA	
LIECHTENSTEIN	
LITHUANIA	
MACEDONIA	
MADAGASCAR	

07

ROOT CANAL
OBTURATION

- Millimeter Marked Gutta Percha Points
- Special Taper Gutta Percha Points
- Millimeter Marked Paper Points
- Special Taper Paper Points
- Temporary Filling Materials
- Diapex Plus (Premixed Calcium Hydroxide Paste with Iodoform)
- Diapaste (Premixed Calcium Hydroxide Paste with Barium Sulfate)
- DiaPen(Cordless Warm Vertical Compaction Device)
- DiaGun(Cordless Backfill Obturation Device)

22

COSMETIC &
RESTORATIVES

- Restorative Micro Hybrid Composite Resin
- Flowable Restorative Composite Resin
- Light Cure Bonding Agent
- Etching Gel
- Light Cure Pit & Fissure Sealant
- Dual Cure Core Build-up Composite
- Self Curing Core Build-up Composite
- Cordless LED Curing Light
- Temporary Dental Cement
- EDTA Root Canal Preparation Cream

35

ENDODONTIC
INSTRUMENTS

- Nickel Titanium Rotary Files
- Flexible K-File
- Lentulo
- Probe File
- Gate Drills
- Peeso Reamers
- Nickel Titanium Hand Files
- Stainless Steel Hand Files
- Stainless Steel Spreader
- Stainless Steel Plugger
- Endo Condenser
- Endo Gauges
- Irrigation Syringes
- Endodontic Irrigation Tips

45

MISCELLANEOUS
ITEMS

- Aluminum Endo Storage Boxes
- Protective Shields & Glasses
- Sterilizing Containers & Trays
- Prophy Brushes
- Prophy Cups
- Bur Blocks & Stands
- Instrument Storage Trays
- Nickel Titanium Arch Wires
- Glass Items
- Dental Wax

ROOT CANAL
OBTURATION


- Millimeter Marked Gutta Percha Points
- Special Taper Gutta Percha Points
- Millimeter Marked Paper Points
- Special Taper Paper Points
- Temporary Filling Materials
- Diapex Plus (Premixed Calcium Hydroxide Paste with Iodoform)
- Diapaste (Premixed Calcium Hydroxide Paste with Barium Sulfate)
- DiaPen(Cordless Warm Vertical Compaction Device)
- DiaGun(Cordless Backfill Obturation Device)

DiaDent’s Millimeter Marked Gutta Percha Points & Paper Points


- All DiaDent Points are Hand Jig Rolled to exacting standards.
- They are uniformly tapered and unsurpassed for consistency in sizing.
- All Points are laser inspected to ensure the tightest tolerances surpassing ISO standards.

Millimeter Marked Gutta Percha Points

- Quick and easy method for depth measurements.
- Points with mm marking at 16mm, 18mm, 19mm, 20mm, 22mm and 24mm.
- Eliminates apex perforations which may lead to bleeding.
- Reaches the apex, completely seal the entire length of the canal.
- Save chair time with premeasured points.
- Most accurate depth measurements.



Millimeter Marked




Computerized Laser Inspection System

- Each point is inspected by the laser beam.
- Using a laser beam, each DiaDent ML.029 point is individually inspected for diameters at D3 and D16.
- ISO and ADA standard requires the use of 1/100mm incremental measurement units, DiaDent uses 1/1,000mm increments.


Millimeter Marked Paper Points

- Quick and easy method for depth measurements.
- Points with mm marking at 18mm, 19mm, 20mm, 22mm and 24mm.
- Eliminates apex perforations which may lead to bleeding.
- Completely drying the entire length of the canal.
- Save chair time with premeasured points.
- Most accurate depth measurements.




Millimeter Marked

DiaDent Quality Assurance



Comparative absorption of Paper Points, by Brand




DIADENT PATENTS

- US Patent No. 6,004,133
- International Patents pending.
- European Patents No. 0964656
- United Kingdom Patents No. 0964656
- Deutsches Patents No. EP 0964656

Hand Jig Rolled Points vs Machine Rolled Points


- DiaDent’s method of hand rolling the points ensures the following:

DiaDent Hand Jig Rolled Gutta Percha Points



- Consistent, Uniform taper
- Smooth, Rounded Bullet shaped tips
- Readily conforms to shape of canal
- Greater likelihood of optimal apical seal

Machine Rolled Gutta Percha Points



- Inconsistent taper
- Uneven, Unpredictably shaped tips
- Overlapping of material of the tips
- Reduced likelihood of optimal apical seal

ML.029

MILLIMETER MARKED + LASER INSPECTED GUTTA PERCHA POINTS



ML.029

- M - Millimeter Marked
- L - Laser Inspected
- .029 - Diameter Tolerance

Available Sizes

Size	P-63	SP-60
XF	ML 102-601	ML 102-S601
FF	ML 102-602	ML 102-S602
MF	ML 102-603	ML 102-S603
F	ML 102-604	ML 102-S604
FM	ML 102-605	ML 102-S605
M	ML 102-606	ML 102-S606
ML	ML 102-607	ML 102-S607
L	ML 102-608	ML 102-S608
XL	ML 102-609	ML 102-S609
XF/M	ML 102-691	ML 102-S691

(Pkg. of 100)

Package Information

• P-43 Vial Package

Contains 6 light weight plastic square vials with color coded caps
Box dimension :
77 × 17 × 50mm(W × H × D)

• P-63 Spill Proof Package

6 compartments with outer paper box
Box dimension:
73 × 17 × 49mm(W × H × D)

• SP-60 Slide Package

one hand slide open lid
Box dimension:
70 × 13 × 45mm(W × H × D)

Size vs Color Chart

15	45	90	XF	ML
20	50	100	FF	L
25	55	110	MF	XL
30	60	120	F	08
35	70	130	FM	10
40	80	140	M	

(Pkg. of 120)

COLOR CODED GUTTA PERCHA POINTS



Available Sizes

Size	P-60
XF	102-601
FF	102-602
MF	102-603
F	102-604
FM	102-605
M	102-606
ML	102-607
L	102-608
XL	102-609
XF/M	102-691

(Pkg. of 100)

Package Information

• P-40 Vial Package

Contains 6 light weight plastic round vials with color coded caps
Box dimension :
79 × 16 × 53mm(W × H × D)

• P-60 Spill Proof Package

6 compartments plastic package
Box dimension :
71 × 13 × 46mm(W × H × D)

Size vs Color Chart

15	45	90	XF	ML
20	50	100	FF	L
25	55	110	MF	XL
30	60	120	F	08
35	70	130	FM	10
40	80	140	M	

(Pkg. of 120)

DIA-PROISO.04/.06 PLUS MILLIMETER MARKED + LASER INSPECTED SPECIAL TAPER GUTTA PERCHA POINTS

- Millimeter Marked Special Taper Gutta Percha Points.
- Designed to be used to complement for **Profile, Hero & K3** files.
- Available in taper .04 & .06



Available Sizes

Size	Dia-ProISO .04 Plus		Dia-ProISO .06 Plus	
	P-63	SP-60	P-63	SP-60
#15	ML 158-603	ML 158-S603	ML 159-603	ML 159-S603
#20	ML 158-604	ML 158-S604	ML 159-604	ML 159-S604
#25	ML 158-605	ML 158-S605	ML 159-605	ML 159-S605
#30	ML 158-606	ML 158-S606	ML 159-606	ML 159-S606
#35	ML 158-607	ML 158-S607	ML 159-607	ML 159-S607
#40	ML 158-608	ML 158-S608	ML 159-608	ML 159-S608
#45	ML 158-609	ML 158-S609	ML 159-609	ML 159-S609
#50	ML 158-610	ML 158-S610	ML 159-610	ML 159-S610
#55	ML 158-611	ML 158-S611	ML 159-611	ML 159-S611
#60	ML 158-612	ML 158-S612	ML 159-612	ML 159-S612
#70	ML 158-613	ML 158-S613	ML 159-613	ML 159-S613
#80	ML 158-614	ML 158-S614	ML 159-614	ML 159-S614
#15/40	ML 158-691	ML 158-S691	ML 159-691	ML 159-S691
#45/80	ML 158-692	ML 158-S692	ML 159-692	ML 159-S692

(Pkg. of 60)

DIA-PROT PLUS



- Designed to be used to complement **Protaper** Files.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#F1	MLP 150-S601
#F2	MLP 150-S602
#F3	MLP 150-S603
#F4	MLP 150-S604
#F5	MLP 150-S605
#F1/F3	MLP 150-S691
#F4/F5	MLP 150-S692
#F1/F5	MLP 150-S693

(Pkg. of 60)

DIA-PRO R



- Designed to be used to complement **RECIPROC®** Files.
- Millimeter Marked Points

Size	ITEM No.
#25	ML164-605
#40	ML164-608
#50	ML164-610
#Asst(25/40/50)	ML164-691

(Pkg. of 60)

DIA-PRO W



- Designed to be used to complement **Wave One** Files.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#Small	ML165-601
#Primary	ML165-602
#Large	ML165-603
#Asst(Small/Primary/Large)	ML165-691

(Pkg. of 60)

DIA-T



- Designed to be used to complement **TF** Files.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
.06 #25	ML163-0605
.06 #30	ML163-0606
.06 #35	ML163-0607
.04 #40	ML163-0408
.06 #25/.04 #40	ML163-1691

(Pkg. of 60)

FEATHERED TIP




- Designed to be used as auxiliary points to the master cone technique
- Extra fine taper
- Hand Jig Rolled

Size	ITEM No.
#XF	126-601
#FF	126-602
#MF	126-603
#F	126-604
#FM	126-605
#M	126-606
#ML	126-607
#L	126-608
#XF/M	126-691

(Pkg. of 100)

DIA-GTP



- Designed to be used to complement **Profile GT Greater Taper** Files.
- Available in .06, .08, .10, .12, taper.

Size	ITEM No.
DIA-GTP .06 #20	112-601
DIA-GTP .08 #20	112-602
DIA-GTP .10 #20	112-603
DIA-GTP .12 #20	112-604
DIA-GTP .06/.12 #20	112-691

(Pkg. of 60)

MMPP

MILLIMETER MARKED PAPER POINTS



- MMPP is marked to provide preliminary guidance in measuring the root canal depth before proceeding clinically with the application of the Gutta Percha Points.
- These points are highly absorbent and color coded to identify ISO tip sizes.

Available Sizes

Size	P-33	P-63	SP-60
XXF	MP 202-301	MP 202-601	MP 202-S601
XF	MP 202-302	MP 202-602	MP 202-S602
F	MP 202-303	MP 202-603	MP 202-S603
M	MP 202-304	MP 202-604	MP 202-S604
C	MP 202-305	MP 202-605	MP 202-S605
XC	MP 202-306	MP 202-606	MP 202-S606
XXF/XC	MP 202-391	MP 202-691	MP 202-S691

(Pkg. of 200)

Size vs Color Chart

8	10		
15	45	90	XXF
20	50	100	XF
25	55	110	F
30	60	120	M
35	70	130	C
40	80	140	XC

Package Information

- **P-33 Cell Package**
Points remain sterile until opened
Box dimension: 78 x 19 x 50mm(W x H x D)

- **P-43 Vial Package**
Contains 6 light weight plastic square vials with color coded caps
Box dimension: 77 x 17 x 50mm(W x H x D)

- **P-63 Spill Proof Package**
6 compartments plastic package
Box dimension: 73 x 17 x 49mm(W x H x D)

- **SP-60 Slide Package**
One hand slide open lid
Box dimension: 70 x 13 x 45mm(W x H x D)

Contents

Size	P-33	P-43	P-63	SP-60
#08-80	200pcs	200pcs	200pcs	200pcs
#90-140	200pcs	120pcs	120pcs	120pcs

COLOR CODED PAPER POINTS

- Made from specially formulated highly absorbent papers.
- Color Coded.
- Hand Jig rolled to ensure Superior absorbency.

Available Sizes

Size	P-60
XXF	202-601
XF	202-602
F	202-603
M	202-604
C	202-605
XC	202-606
XXF/XC	202-691

(Pkg. of 200)

Size vs Color Chart

8	10		
15	45	90	XXF
20	50	100	XF
25	55	110	F
30	60	120	M
35	70	130	C
40	80	140	XC

Contents

Size	P-40	P-60
#08-80	200pcs	200pcs
#90-140	120pcs	120pcs

Size	P-40	P-60
#8	201-401	201-601
#10	201-402	201-602
#15	201-403	201-603
#20	201-404	201-604
#25	201-405	201-605
#30	201-406	201-606
#35	201-407	201-607
#40	201-408	201-608
#45	201-409	201-609
#50	201-410	201-610
#55	201-411	201-611
#60	201-412	201-612
#70	201-413	201-613
#80	201-414	201-614
#90	201-415	201-615
#100	201-416	201-616
#110	201-417	201-617
#120	201-418	201-618
#130	201-419	201-619
#140	201-420	201-620
#15/40	201-491	201-691
#45/80	201-492	201-692
#90/140	201-493	201-693

DIA-PROISO .04/.06 PLUS

MILLIMETER MARKED SPECIAL TAPER PAPER POINTS

- Millimeter marked Special Taper Absorbency Paper Points.
- Designed to be used for drying and cleaning the canals prepared with **Profile, Hero & K3 files**.
- Available in .04 & .06 taper.




Available Sizes

Size	Dia-Pro Iso .04 Plus		Dia-Pro Iso .06 Plus	
	P-63	SP-60	P-63	SP-60
#15	MP 258-603	MP 258-S603	MP 259-603	MP 259-S603
#20	MP 258-604	MP 258-S604	MP 259-604	MP 259-S604
#25	MP 258-605	MP 258-S605	MP 259-605	MP 259-S605
#30	MP 258-606	MP 258-S606	MP 259-606	MP 259-S606
#35	MP 258-607	MP 258-S607	MP 259-607	MP 259-S607
#40	MP 258-608	MP 258-S608	MP 259-608	MP 259-S608
#45	MP 258-609	MP 258-S609	MP 259-609	MP 259-S609
#50	MP 258-610	MP 258-S610	MP 259-610	MP 259-S610
#55	MP 258-611	MP 258-S611	MP 259-611	MP 259-S611
#60	MP 258-612	MP 258-S612	MP 259-612	MP 259-S612
#70	MP 258-613	MP 258-S613	MP 259-613	MP 259-S613
#80	MP 258-614	MP 258-S614	MP 259-614	MP 259-S614
#15/40	MP 258-691	MP 258-S691	MP 259-691	MP 259-S691
#45/80	MP 258-692	MP 258-S692	MP 259-692	MP 259-S692

(Pkg. of 100)

DIA-PROISO .04/.06 PLUS
IN THE CELL PACKAGE



- Points remain sterile until opened

Cell pkg	Dia-Pro Iso .04 Plus	Dia-Pro Iso .06 Plus
#15	MP 258-703	MP 259-703
#20	MP 258-704	MP 259-704
#25	MP 258-705	MP 259-705
#30	MP 258-706	MP 259-706
#35	MP 258-707	MP 259-707
#40	MP 258-708	MP 259-708
#15/40	MP 258-791	MP 259-791

(Pkg. of 120)

DIA-PROT




- Designed to be used to complement **Protaper Files**.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#F1	MP250-S601
#F2	MP250-S602
#F3	MP250-S603
#F4	MP250-S604
#F5	MP250-S605
#F1/F3	MP 250-S691
#F4/F5	MP 250-S692
#F1/F5	MP 250-S693

(Pkg. of 100)

DIA-PRO R




- Designed to be used to complement **RECIPROC® Files**.
- Millimeter Marked Points

Size	ITEM No.
#25	MP264-605
#40	MP264-608
#50	MP264-610
#Asst(25/40/50)	MP264-691

(Pkg. of 100)

DIA-PRO W




- Designed to be used to complement **Wave One Files**.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
#Small	MP265-601
#Primary	MP265-602
#Large	MP265-603
#Asst(Small/Primary/Large)	MP265-691

(Pkg. of 100)

DIA-T




- Designed to be used to complement **TF Files**.
- Millimeter Marked Points
- One hand slide open lid

Size	ITEM No.
.06 #25	MP263-0605
.06 #30	MP263-0606
.06 #35	MP263-0607
.04 #40	MP263-0408
.06 #25/.04 #40	MP263-1691

(Pkg. of 100)

FEATHERED TIP



- Designed to be used as auxiliary points to the master cone technique.
- Extra fine taper.
- Hand Jig Rolled.

Size	ITEM No.
#XXF	204-601
#XF	204-602
#F	204-603
#M	204-604
#C	204-605
#XC	204-606
#XXF/XC	204-691

(Pkg. of 200)

ENDO ORGANIZERS

- Packed in plastic Endo Organizer package with 12 compartments.
- The ideal way to organize Endodontic Points.
- DiaDent Endo Organizers can be ordered empty or a combination of Gutta Percha and Paper Points



mmGP1



mmPP2



mmCOMBI1



mmGP2(6%)



mmPP1(4%)



mmCOMBI1(4%)

ML.029 & MMPP

- Millimeter marked Gutta Percha and Paper Points.
- All DiaDent points are Hand Jig Rolled.

Size	Description	ITEM No.
mmEG1	Empty box	EMG 301-100
mmGP1	240points per box(size15/40)	EMG 301-101
mmGP2	240points per box(size45/80)	EMG 301-102
mmGP3	240points per box(size15/80)	EMG 301-103
mmPP1	400points per box(size15/40)	EMP 301-201
mmPP2	400points per box(size45/80)	EMP 301-202
mmPP3	400points per box(size15/80)	EMP 301-203
mmCOMBI1	GP 120points/PP 200points per box(size15/40)	EMC 301-301
mmCOMBI2	GP 120points/PP 200points per box(size45/80)	EMC 301-302
mmCOMBI3	GP 240points/PP 400points per box(size15/40)	EMC 301-303
mmCOMBI4	GP 240points/PP 400points per box(size45/80)	EMC 301-304

.04 & .06 Tapered ML.029 & MMPP

- Millimeter marked special taper Gutta Percha and Paper Points.
- Designed to be used to complement for **Profile, Hero & K3 files**.

Size	Description	ITEM No.
mmGP 1(4%)	240points GP DIA-PROISO-.04 PLUS(size 15/80)	MSG 301-101
mmGP 2(6%)	240points GP DIA-PROISO-.06 PLUS(size 15/80)	MSG 301-102
mmGP 3(4% & 6%)	240points GP DIA-PROISO-.04 & .06 PLUS (size 15/40)	MSG 301-103
mmPP 1(4%)	240points PP DIA-PROISO-.04 PLUS(size 15/80)	MSP 301-201
mmPP 2(6%)	240points PP DIA-PROISO-.06 PLUS(size 15/80)	MSP 301-202
mmPP 3(4% & 6%)	240points PP DIA-PROISO-.04 & .06 PLUS (size 15/40)	MSP 301-203
mmCOMBI 1 (4%)	240points GP & PP DIA-PROISO-.04 PLUS(size15/40)	MSC 301-301
mmCOMBI 2 (6%)	240points GP & PP DIA-PROISO-.06 PLUS(size15/40)	MSC 301-302

DIA-GP GAUGE

- ITEM #781-001**
 - Designed to size a perfect master cone from any Gutta Percha Points.
 - Drop it into the apical size and cut off with a blade.
 - Autoclavable
 - Dimensions: 92 × 3 × 20m/m(W × H × D)



PACKAGE INFORMATION

- P-33 Cell Package**
 - Points remain sterile until opened
 - Box dimension: 78 × 19 × 50mm(W × H × D)
- P-43 Vial Package**
 - Contains 6 light weight plastic square vials with color coded caps
 - Box dimension: 77 × 17 × 50mm(W × H × D)




- P-63 Spill Proof Package**
 - 6 compartments with outer paper box
 - Box dimension: 73 × 17 × 49mm(W × H × D)
- SP-60 Slide Package**
 - One hand slide open lid
 - Cellophane shrink-wrapped package
 - Box dimension: 70 × 13 × 45 mm (W × H × D)



TEMPORARY STOPPING

- For temporary filling and thermal diagnostics tool.
- A warmed stick with pinpoint thermal sensitivity.
- High quality Gutta Percha.
- Easy to heat.
- Available in White

	White
120g	180-201
250g	180-301




DIATEMP™ A VISIBLE, LIGHT CURED TEMPORARY FILLING MATERIAL

Light Cured material with nano-sized silver for temporary/intermediate restorations of inlay and onlay, for relining prefabricated crowns and bridges.



- Benefits**
- Quick and simple removal: Can be removed in one piece
 - High wear resistance, mechanical strength
 - Ensures a fast set
 - Excellent adhesion
 - Easy application: Convenient syringe dispensed material does not required mixing
 - Blue Dia-Temp is ideal for Posterior teeth
 - Yellow Dia-Temp is ideal for Anterior teeth

- Packages**
- 3 syringes of 3g
 - ITEM #2003-1101 Blue
 - ITEM #2003-1201 Yellow



E-TEMP HYDRAULIC TEMPORARY RESTORATIVE MATERIAL

- Temporary filling material for standard cavities and post-endodontic treatments.



Benefits


- Temporary Filling Material for Standard Cavities.
- Easily setting initiated by saliva.
- Excellent adaptability.
- Easy to remove.
- Adheres to dentine, expands slightly during the setting and hardens rapidly to ensure hermetical seal.
- Improved the mechanical properties for the quick temporary filling.

Directions


- Prepare and dry the cavity with a small cotton pellet.
- Transfer E-TEMP to the cavity using a suitable instrument and form the contour. In the case of occlusal Class I, the patient bites down on the surface, and then form the contour.
- After the patient close the mouth, E-TEMP will harden within 25 minutes. Recommend the patient to avoid the strong pressure or abrasion for at least 2 hours.

Package

- 30g ITEM #2004-1001



- This test is measured as the same condition at 37 degrees and 95% humidity in the laboratory(R&D DiaDent 2013)





- This test is measured as ISO 3107 regulations after the same setting time in the laboratory(R&D DiaDent 2013)

DIAPEX® PLUS

CALCIUM HYDROXIDE PASTE WITH IODOFORM

- Premixed Calcium Hydroxide + Iodoform paste for root canal treatment.
- Temporary root canal filling material after a pulpectomy.
- Ideal for the treatment of infected root canal and for vital pulpotomies in deciduous teeth.



Indications

- Intracanal medicament
- Apexification
- Exudation Control
- Periapical Lesions
- Root Resorption
- Temporary root filling
- Perforations
- Under developed Pulpless Teeth

Benefits

- Quiets "Hot Abscess"
- Disinfects canal
- Promotes Apexification
- Treats Traumatic Injuries
- Radiopaque
- Quick and simple application
- Antibacterial and Bacterostatic


Packages

- **ITEM #1001-501 Regular Kit Type A**
- 2g of paste and 20 disposable tips with cap
- **ITEM #1001-502 Regular Kit Type B**
- 2g of paste and 10 disposable tips with cap
- **ITEM #1001-503 Regular Kit Type C**
- 2g of paste and 4 disposable tips with cap
- **ITEM #1001-504 Refill Kit**
- 2g of paste
- **ITEM #1001-505 Intro Kit**
- 0.5g of paste and 4 disposable tips
- **ITEM #1001-316 Disposable Syringe Tips**
- 20 disposable tips in a **blister package**

- This Item is not available in Europe

DIAPASTE™

CALCIUM HYDROXIDE PASTE WITH BARIUM SULFATE



Characteristics

- Premixed Calcium Hydroxide + Barium Sulfate paste for root canal treatment.
- Temporary root filling.
- Root canal filling for primary teeth.
- Use after pulp extirpation treatment.
- Filling on the pulpotomy.
- Temporary Pulp Capping.

Benefits

- Water Soluble: Easy to clean and remove, reduce operative time
- Easy and quick access
- Radiopaque
- Antibacterial

Packages

- **ITEM #1001- 401 Regular Kit Type A**
- 2g of paste and 20 disposable tips with cap
- **ITEM #1001- 402 Regular Kit Type B**
- 2g of paste and 10 disposable tips with cap
- **ITEM #1001- 403 Regular Kit Type C**
- 2g of paste and 4 disposable tips with cap
- **ITEM #1001- 404 Refill Kit**
- 2g of paste
- **ITEM #1001- 405 Intro Kit**
- 0.5g of paste and 4 disposable tips
- **ITEM #1001- 315 Blue Disposable Syringe Tips**
- 20 disposable tips in a **blister package**

DIAPEN CORDLESS WARM VERTICAL COMPACTION DEVICE

- DiaDent's Cordless Warm Vertical Compaction Device
- DiaPen has the optimum functionality and quality and is fully proven effective for the downpack obturation system. DiaPen has the following techniques and features:

Techniques

- **Softening** the root canal filling material(Gutta Percha Points) during the root canal treatment
- **Spreading** the root canal filling material(Gutta Percha Points) laterally and vertically
- **Cutting** the root canal filling material(Gutta Percha Points)
- **Compacting** the root canal filling material(Gutta Percha Points) using with the heat

Features

- Excellent result with warm and vertical compaction for all canals including lateral canals. DiaPen offers reliable results with tightly sealed root canals
- Cordless, Ergonomic and Lightweight design make it easy for handling, cleaning and storage
- Quick heating tip reaches high level of temperature within one (1) second to save treatment time
- Color-coded pen tips are available in various sizes (#XF, F, FM, M, ML) for easy identification
- Three temperature settings (High, Medium, Low)
- Simple operation


Package

- **ITEM #4002-1001 Regular Kit**
- DiaPen Handpiece & Charger
- Accessories: Silicone Cover(2ea), Power Adapter & Cord, Pen Tips #F (included Endo stops, 1ea), #FM(included Endo Stops, 1ea)

Accessories

- **ITEM #4002-1010 Handpiece**
- **ITEM #4002-1020 Charger**
- **ITEM #4002-1030 Power Adapter & Cord**
- **ITEM #4002-1040 Silicone Cover(2ea)**
- **ITEM #4002-1050 Battery**
- **ITEM #4002-1060 Pen Tip #XF(1ea)**
- **ITEM #4002-1061 Pen Tip #F(1ea)**
- **ITEM #4002-1062 Pen Tip #FM(1ea)**
- **ITEM #4002-1063 Pen Tip #M(1ea)**
- **ITEM #4002-1064 Pen Tip #ML(1ea)**
- **ITEM #4002-1065 Pen Tips #F/FM(1ea each)**

Color-coded pen tips: #XF, F, FM, M, ML



DIAGUN CORDLESS BACKFILL OBTURATION DEVICE

- DiaGun is Cordless Obturation System Device, DiaGun offers fast and easy obturation delivery. With this effective obturation system and technique, optimum results are achieved. The high-performance DiaGun device has the following features:

Features

- Cordless and Compact Design allows complete freedom of movement
- Reliable Gutta Percha Filling to the desired level without delay
- Fast and continuous flow of filling tightly seals the canal
- Variable Temperature Settings (160°C, 180°C, 200°C): precise control of obturation flow
- Quick Heating System: temperature reaches 200°C within 25 seconds
- Ergonomically designed 360° Swivel Needle offers improved access
- Thinner needle can easily fill the narrow root canal

Package

- **ITEM #4001-1001 Regular Kit**
- DiaGun Handpiece & Charger
- Accessories: Safety Cap(2ea), Plunger, Plunger O-Rings(4ea), Power Adapter & Cord, Disposable Gun Tips(23G,4ea), Multipurpose wrench, Scrubbing Brush, Battery, Gutta Percha Obturator(100ea)

Accessories

- **ITEM #4001-1010 Handpiece**
- **ITEM #4001-1020 Charger**
- **ITEM #4001-1030 Safety Cap(2ea)**
- **ITEM #4001-1040 Plunger**
- **ITEM #4001-1050 Plunger O-Rings(4ea)**
- **ITEM #4001-1060 Power Adapter & Cord**
- **ITEM #4001-1070 Disposable Gun Tips(23G,4ea)**
- **ITEM #4001-1071 Disposable Gun Tips(25G,4ea)**
- **ITEM #4001-1080 Multipurpose Wrench**
- **ITEM #4001-1090 Scrubbing Brush**
- **ITEM #4001-1100 Battery**
- **ITEM #181-101 Gutta Percha Obturator(Medium type, 100ea)**



GUTTA PERCHA OBTURATOR


- **ITEM #181-100 Soft type**
- **ITEM #181-101 Medium type**
- **ITEM #181-102 Hard type**

- Designed to use with DiaGun the obturation gun.
- Specially formulated Gutta Percha with excellent flowing characteristics.
- Easy to heat
- Pkg. of 100




ROOT CANAL OBTURATION TECHNICAL SEQUENCE

DOWNPACK OBTURATING WITH DIAPEN




ROOT CANAL OBTURATION TECHNICAL SEQUENCE

BACKFILL OBTURATING WITH DIAGUN



Radiographs



Perfect and tightly sealed root canals after the root canal warm obturation technique with DiaPen & DiaGun.

COSMETIC & RESTORATIVES

- Restorative Micro Hybrid Composite Resin
- Flowable Restorative Composite Resin
- Light Cure Bonding Agent
- Etching Gel
- Light Cure Pit & Fissure Sealant
- Dual Cure Core Build-up Composite
- Self Curing Core Build-up Composite
- Cordless LED Curing Light
- Temporary Dental Cement
- EDTA Root Canal Preparation Cream

DIAFIL™ A LIGHT CURE, ESTHETIC RESTORATIVE MICRO HYBRID COMPOSITE RESIN

- Micro hybrid composite resin for use in Anterior and Posterior restorations



Start Kit

Benefits

- Unique fillers offer minimum polymerization shrinkage which helps reduce stress and sensitivity on the tooth potentially
- Has excellent fracture toughness, high tensile, compressive strengths for long lasting results
- Very easy handling of acceptable viscosity and good marginal adaptation
- Remarkable color stability and superior quality of esthetics
- High Level of radiopacity

Packages

- **Start Kit**
 - 5 composite resin syringes (4g)
 - 1 Bottle bonding agent(5ml)
 - 1 Etching gel syringe(3ml)
 - Accessories: 50 brush tips, 1 brush handle, mixing well, mixing pad, 10 disposable needle tips, shade sheet

Set	Composite Resin Shades	ITEM No.
START KIT A SET	A1, A2, A3, A3.5, B2	2001-1101
START KIT B SET	A2, A3, A3.5, B2, B3	2001-1102
START KIT C SET	A1, A2, A3, A3.5, C2	2001-1103



Refill Package

- **Refill Package**
 - 1 syringe x 4g

A1	2001-1201	B2	2001-1207
A2	2001-1202	B3	2001-1208
A3	2001-1203	C2	2001-1210
A3.5	2001-1204	C3	2001-1211
A4	2001-1205	A20	2001-1212
B1	2001-1206	A30	2001-1213





Capsule Type

- **Capsule Type**
 - 0.25g x 20 capsules



A1	2001-1405
A2	2001-1401
A3	2001-1402
A3.5	2001-1403
B2	2001-1404

PHYSICAL PROPERTIES


Flexural Strength




Compressive Strength



Depth of Cure



Abrasive, Wear Resistance




• POSTERIOR RESTORATION




• ANTERIOR RESTORATION




• HIGH RADIOPACITY




TECHNIQUE GUIDES




1. Prepare the tooth : Isolate the tooth with rubber dam and clean the tooth




2. Etch : Apply DiaEtch(etchant) to enamel and dentin and wait 15 seconds. Rinse and blot the excess water, leaving tooth moist




3. Bond : Apply 3 consecutive coats of DiaPlus (bonding agent) to enamel and dentin and wait 15 seconds. Dry gently for 3 seconds and light cure for 20 seconds



4. Place composite resin : Place DIAFIL™(composite resin) in layers less than 2.0mm



5. Light Cure : Cure each increment for 20 seconds



6. Finish and Polish

DIAFIL™ FLOW LIGHT CURED FLOWABLE RESTORATIVE COMPOSITE RESIN

Flowable composite with very easy handling and high strength



Benefits

- Has a high level of fluoride release with no slump, run which helps a fast techniques and handling
- Enables to be used in various restorative procedures
 - : Class III, V and IV direct restorations
 - A liner for Class I, II
 - Pit & Fissure sealant where is required
 - Well suited for minimally invasive dentistry
- Excellent tensile strength and low shrinkage
- High radiopaque materials

Packages



- **Economic Package**
 - Flowable Resin 2g x 4 syringes (Available shade A1, A2, A3, A3.5, B2, B3, C2, A20, A30)
 - 40 disposable tips(black)
 - Shade Sheet
- **Refill Package**
 - Flowable Resin 2g x 1 syringe (Available shade A1, A2, A3, A3.5, B2, B3, C2, A20, A30)
 - 10 disposable tips(black)

SHADE	Economic Pkg.	Refill Pkg.
A1	2001-4101	2001-4201
A2	2001-4102	2001-4202
A3	2001-4103	2001-4203
A3.5	2001-4104	2001-4204
B2	2001-4105	2001-4205
B3	2001-4106	2001-4206
C2	2001-4107	2001-4207
A20	2001-4108	2001-4208
A30	2001-4109	2001-4209

PHYSICAL PROPERTIES



Flexural Strength

The flexural strength is measured by pressing a stick-typed experimental material. Flexural strength is the value obtained when the sample breaks, fixing both ends of sample and pressing in the middle of it. This test combines the forces found in compression and tension. DIAFIL™ FLOW shows statistically higher than other 2 brands.





Compressive Strength

Compressive strength is particularly important because it is related with masticating forces. Rods are made of the material and simultaneous forces are applied to the opposite ends of the sample length vertically. The sample failure is a result of shear and tensile forces. DIAFIL™ FLOW is statistically higher than other brands.



Flexural Modulus

Flexural modulus is a method of defining a material's stiffness. A low modulus indicates high flexible material. DIAFIL™ FLOW shows a similar flexural elastic rate with other flowable composite resins.



• RESTORATION WITH DIAFIL™ FLOW



DIAPLUS™ SINGLE COMPONENT BONDING AGENT


Recommended for all direct composite restorations, light-cured bonding procedures



Benefits

- High dentin, enamel bond strength performance
- Available in all classes of direct composite restoration
- Indirect restorations of Porcelain, Amalgam or Composite Repair after light curing with a minimal film thickness
- Ethanol-based bonding agent enables to maintain its exceptional bond strength

Shear bond strength(MPa)



Packages

- **ITEM #2001-2101**
- one-bottle bonding agent (1 bottle x 5ml)

DISPOSABLE BRUSHES

- Flexible adhesive tip allows restorations and small objects of all kind to be picked up, held, carried and placed quickly and easily
- These brushes are included DiaPlus G7's Regular Package

Size	ITEM NO.
Small	2001-6301
Medium	2001-6302
Large	2001-6303

(Pkg. of 100)



DIAPLUS™ G7 SINGLE COMPONENT SELF-ETCHING LIGHT-CURED BONDING AGENT

DiaPlus™ G7, the innovative 7th Generation bonding agent and offers superior and consistent bonding results on dentin and enamel with one bottle, one coat self-etching bonding agent system.



Benefits

- Has **Etch, Desensitize, Prime and Bond** all with One-Bottle, One-Coat bonding system and no separate etching
- **Fast 3 Steps** (less than 35 seconds) reduces the risk of contamination
- Offer a **tight and long-lasting bond** between light-cured composite and both enamel and dentin by a highly advanced formulation whether wet or dry on dentin and enamel
- **Strong and tight bond** to dentin and enamel is offered by 4-MET monomer for adhesive to dentin and Phosphoric Acid Ester Monomer for superior bonding strength to enamel
- Consistent and a **long term bond strength** with the unique Nano-Filler
- Provide a very low risk of post-operative sensitivity
- **Remarkable thin bonding layer** less than 10^{um} offers a frosty surface for much easier application of the first composite layer


Packages

Regular Package


- **ITEM #2001-2701**
- One bottle bonding agent (1 bottle x 5ml)
- 1 Mixing well
- 50 disposable brushes

Refill Package

- **ITEM #2001-2702**
- One bottle bonding agent (1 bottle x 5ml)




Thin Bonding Layer



Superior Bond Strength

- One coat, simple and fast bonding application



Less than 35 seconds

Technique Guide

DIAETCH™ 37%

A 37% ETCHING GEL PHOSPHORIC ACID

A phosphoric acid is used in etch enamel and/or dentin

Benefits

- Phosphoric acid etching of enamel and dentin with ideal thickness
- Quick and effective application to the desired place
- Enables to rinse completely and easily from the tooth surface
- High biocompatibility

Packages

Economic Package Type A

- **ITEM #2001-3101** - 5 Syringes × 5ml and 20 disposable tips

Economic Package Type B

- **ITEM #2001-3102** - 5 syringes × 3ml and 20 disposable tips

Regular Package

- **ITEM #2001-3103** - 2 syringes × 3ml and 10 disposable tips

Refill Package Type A

- **ITEM #2001-3202** - 1syringe × 5ml and 5 disposable tips



Before After
(Microphotographs of premolar, X1,000)

DIASEAL™

A LIGHT CURE PIT & FISSURE SEALANT

Designed for sealing enamel pits and fissures of teeth to prevent caries

Benefits

- UDMA series monomer and Nano-filler offers high flowability, excellent wear resistance and strength
- Optimal viscosity to flow easily and helps for easy application
- Ensures hygienic and easy delivery and outstanding versatility
- High visual verification during application

Packages

Regular Package

- **ITEM #2001-5101**
 - Sealant 2 syringes × 1.2ml
 - Etching 2 syringes × 1.2ml
 - 10 disposable needle tips for Dia-Seal
 - 10 disposable needle tips for Dia-Etch

Refill Package

- **ITEM #2001-5201**
 - Sealant 1 syringe × 1.2ml
 - 5 disposable needle tips

Properties

Mean value and standard deviation of micro-leakage score

Materials	Micro-leakage score (standard deviation)
DiaSeal™	1.0000(1.1094)
U	1.7333(1.4829)
M	0.9167(1.2795)

Thermocycling 500 times, 1% Methylene blue solution



DIAFIL CORE™ AUTOMIX

DUAL CURED CORE BUILD-UP COMPOSITE

NEW

Dual-curing(Light and self-curing) core build-up and radiopaque two-component composite materials in Automix delivery system

Benefits

- Comfortable dispenser: Automix syringes deliver a correct and consistent mix on every uses
- Easy and precise application with Automix and Intraoral tips
- Dual curing: Light cure for 40 seconds or Chemical cure within 4 minutes
- High compressive strength composite resin and long-term durability
- Cuts and wears like dentin
- Superior bond strength with DiaPlus
- Available 2 colors : Natural color(A3 shade) for the optimal esthetic results and Blue color

Packages

- **ITEM#2001-8101 Natural Color(A3 shade)**
- **ITEM#2001-8102 Blue Color**
 - 2 syringes of 9g
 - 10 Intraoral tips
 - 10 Automix tips
- **ITEM #2001-8109**
 - 20 Automix tips and 20 Intraoral tips



DIAFIL CORE™

SELF-CURING CORE BUILD-UP COMPOSITE

Self-curing and long-term fluoride releasing composite materials consist of two-components, base plus catalyst for core build-up system

Benefits

- Provides good handling properties for easier use
- Offers more faster self cure set time
- Smooth and creamy materials help fill it into needle tips very easily
- High level of radiopacity

Packages

- **ITEM #2001-7101**
 - 10g Base + 10g Catalyst
 - 15 Spatulas
 - 50 Mixing pads
 - 10 Needle tips



D-LUX

CORDLESS LED CURING LIGHT



D-Lux is a new advanced cordless LED curing light. Made from the highest quality aluminum, D-Lux dissipates heat quickly to prevent overheating problem so that it can produce the maximum light intensity safely and consistently.

Features


- **OLED panel:** unique and advanced one-button control
- A high intensity of **1,600 mW/cm²** offers extremely short curing times of no more than **5 seconds**
- Cordless, highly portable and lightweight (only 155g)
- Unique ergonomic design: provides comfortable grip and easy-to-read OLED panel
- Durable one-piece aluminum design is robust and easy to disinfect
- High-performance lithium-ion battery technology:
 - A new fully charged battery provides 160 x 10-second cures
 - OLED panel conveniently indicates battery life and a countdown timer for the next needed charge
 - Battery can be fully recharged in 90 minutes
- **5 versatile curing modes** for every indications
- Easy to change the battery pack

Product Features


Simple & Intuitive OLED Button




High Intensity Output




Fast Charging System (90min)




5 Curing Modes




5 Max:
Maximum light intensity output for 5 seconds




20 Pulse:
Pulses every 0.15 seconds for 20 seconds at maximum light intensity



10 Max:
Maximum light intensity output for 10 seconds



20 Soft:
An initial 10% output of maximum light intensity that gradually increases to full power in 10 seconds, followed by 10 seconds output at full power



20 Mid:
70% output of maximum light intensity for 20 seconds



Package

- **ITEM #4004-1001 Regular Kit**
 - D-Lux Handpiece & Charger
 - Accessories: Light Guide(11 × 11mm), Plastic Protector, Power Supply and Cord

Device Part Codes

- **ITEM #4004-1110 Handpiece**
- **ITEM #4004-1120 Charger**
- **ITEM #4004-1130 Light Guide(11 × 11mm)**
- **ITEM #4004-1150 Battery Pack**
- **ITEM #4004-1160 Plastic Protector**
- **ITEM #4004-1170 Power Supply (Adapter)**
- **ITEM #4004-1171 Power Cord**

DiaCem™ TempNE
NON EUGENOL TEMPORARY CEMENT
NEW

- Non-Eugenol Temporary cement for trial cementing restorations or temporary crowns and bridges



- Features**
 - Easy and precise direct dosing
 - No hand mixing and Stable mix of base and catalyst viscosity with the automix tip
 - High bond strength to prevent leakage and sensitivity to cold and heat
 - Easy removal for the patient comfort
 - Low film thickness with high compressive and flexural strength
 - Fast working time(3 min) and setting time(5 min in the mouth)

- Packages**
 - ITEM#2004-2101 2 Syringes of 10g, 10 Automix tips
 - ITEM#2001-6112 20 Automix tips

DIAPREP™ Plus
EDTA ROOT CANAL PREPARATION CREAM

EDTA Cream with Urea Peroxide, used as an irrigation agent for chemo-mechanical preparation of the root canal.



- Characteristics**
 - Debridement in the root canal
 - Dissolution of pulp remnants
 - Removal of smear layer
 - Facile movements of canal instruments
 - Microbicidal treatment within the canal
 - Clearance of dentin tubules and apical area

- Benefits**
 - EDTA removes smear layer, softens dentin, and facilitates the removal of calcific obstruction.
 - Urea Peroxide reacts with Sodium Hypochlorite to produce oxygen bubbles that flush out debris within the canal.
 - The lubricant medium facilitates instrument movement within the canal to eliminate instrument binding or breakage.
 - Unclogged dentin tubules and apical area permit better penetration of root canal medicaments and fillings.
 - Provides rapid cleansing and enlarging of canals due to EDTA, which softens walls and smooth calcium deposits.
 - Easier placement into canals than liquids.

- Packages**
 - ITEM # 2002-1100 - 2 syringes of 6g, 10 disposable tips and 1 protective cap
 - ITEM # 2001-6111 - 20 Disposable tips in a polybag

ENDODONTIC
INSTRUMENTS

- Nickel Titanium Rotary Files
- Flexible K-File
- Lentulo
- Probe File
- Gate Drills
- Peeso Reamers
- Nickel Titanium Hand Files
- Stainless Steel Hand Files
- Stainless Steel Spreader
- Stainless Steel Plugger
- Endo Condenser
- Endo Gauges
- Irrigation Syringes
- Endodontic Irrigation Tips

HELIX™

NICKEL TITANIUM ROTARY FILES



Benefits



- Made of Superior Nickel Titanium Alloy for Engine Driven use.
- Specially designed for the crown-down technique.
- Offers a smooth and progressive enlargement of the canal.
- Helix NiTi Rotary Files help to ease around curved canals and make a tapered preparation.
- Helix NiTi Rotary Files use ISO Color identification and come in 21mm and 25mm lengths.
- With amazing flexibility, easily move through canals, especially severely curved canals.
- High resistance to cyclic fatigue and file separation.
- Offers outstanding flexibility, cutting ability and reduce user fatigue.
- Ergonomic handle provides comfort and excellent tactile feedback.
- Recommended Speed: 150~350r/min

Direction

- ① Place #10 K-File to as far as it goes in the canal or to approximately Coronal 2/3rd. The file is preferably coated with Dia-Prep Plus (EDTA gel).
- ② Prepare Glide path with #15 K-file to the Level previously explored by #10 K- file.
- ③ Shape with Helix NiTi Rotary files #40 (Black ring) with 0.06 taper (2 lines) to Coronal 1/3rd to Coronal 2/3rd.
- ④ Additionally shape with Helix NiTi Rotary files #30 (Blue Ring) 0.06 taper to Coronal 2/3rd.
- ⑤ Explore with #8 K-File, #10 K-File to the Working Length determined with Apex Locator(NaviRoot Plus) and X-ray.
- ⑥ Prepare Glide Path with #15 K-file to the Working Length.
- ⑦ Shape with Helix NiTi Rotary files #20 (Yellow ring) with 0.04 taper (1 line) to the Working Length.
- ⑧ Determine to what size Apical Shaping is needed.
- ⑨ Shape with Helix Rotary NiTi files #20, #25 or #30 with 0.04 or 0.06 Taper for final Apical Shaping.
- ⑩ Frequent irrigation with Sodium Hypochlorite solution and Lubrication with Dia-Prep Plus (EDTA Gel) throughout the Procedure is recommended.
- ⑪ Proceed to fit appropriate size of DiaDent Gutta Percha Points cones (Dia-ProIso Plus).

Size	21MM		25MM	
	.04	.06	.04	.06
#20	NT505-104	NT506-104	NT505-204	NT506-204
#25	NT505-105	NT506-105	NT505-205	NT506-205
#30	NT505-106	NT506-106	NT505-206	NT506-206
#35	NT505-107	NT506-107	NT505-207	NT506-207
#40	NT505-108	NT506-108	NT505-208	NT506-208
#20/35	NT505-193	NT506-193	NT505-293	NT506-293

(Pkg. of 6)





HELIX ORGANIZER



- **ITEM #806-103 Red**
- **ITEM #806-104 Blue**
 - Ideal for Organizing Helix Files
 - Holds 10 Helix Files
 - Removable insert for Cleaning
 - Dimensions: 50 x 50 x 45mm(W x H x D)

DIA-FLEX FILE

FLEXIBLE K-FILE



- Made of high-quality Stainless Steel for creating a more fracture-resistant shape and edge
- Enables a complete preparation of even severely curved and narrow canals
- Offers better Cutting Power and Flexibility to shape any root canals easily, safely and very fast
- Available in Sizes 15 through 40 and 21mm, 25mm and 31mm lengths

Size	21MM	25MM	31MM
#15	505-103	505-203	505-303
#20	505-104	505-204	505-304
#25	505-105	505-205	505-305
#30	505-106	505-206	505-306
#35	505-107	505-207	505-307
#40	505-108	505-208	505-308
#15/40	505-191	505-291	505-391

(Pkg. of 6)

DIA-SPIRAL FILLER



- Lentulo Paste filler used to evenly distribute root canal sealer and cement in the root canal system
- Conical shape made from grinded wire and made of high grade stainless steel with gradual spirals
- Rather stiff without the safety spiral
- Available in Sizes 25 to 40 and 21mm, 25mm lengths

Size	21MM	25MM
#25	507-105	507-205
#30	507-106	507-206
#35	507-107	507-207
#40	507-108	507-208
#25/40	507-191	507-291

(Pkg. of 4)

DIA-PROBE FILE



- Canal Probe
- The minimal taper helps to perform intermediate steps in the calcified or curved canals
- Available in Size #12 and 21mm, 25mm lengths

Size	21MM	25MM
#12	506-101	506-201

(Pkg. of 6)

GATE DRILLS

- Used for time-saving preparation of the upper and middle part of non-curved canal or the straight part of curved canal effectively and creating the excellent access for Reamers and Files.
- Enlarge the straight coronal portion of the canal with Gate Drills and then step down with Endodontic Files using the biggest size firstly as the Crown Down Technique.
- Superior Stainless Steel with the improved cutting capacity.
- Circle ring is coded on the shank.
- Recommended Speed: 800~1200r/min

Size	28mm	32mm
#1	531-101	531-201
#2	531-102	531-202
#3	531-103	531-203
#4	531-104	531-204
#5	531-105	531-205
#6	531-106	531-206
#Asst1/6	531-191	531-291

(Pkg. of 6)



PEESO REAMERS

- Used for the preparation of coronal third and picking up the fillings inside the root canal.
- Often used for the post path.
- Excellent stainless steel with a sharper blade to provide the rapid enlargement.
- Circle ring is coded on the shank.
- Recommended speed: 800~1200r/min

Size	28mm	32mm
#1	541-101	541-201
#2	541-102	541-202
#3	541-103	541-203
#4	541-104	541-204
#5	541-105	541-205
#6	541-106	541-206
#Asst1/6	541-191	541-291

(Pkg. of 6)



ENDODONTIC NICKEL TITANIUM HAND FILES/PLASTIC PACKAGE

Increased Flexibility

: NiTi Shaft is made of premium Nickel-Titanium with 5 times the flexibility compared to Stainless Steel Files.

Strength

: NiTi has 10 times the stress resistance compared to Stainless Steel Files.

Convenience

: Produces smooth and progressive enlargement of root canal.
NiTi follow existing canal path and has memory so there is no need to Pre-Bend the files.

Efficiency

: Color Coded handles and pre-inserted silicone stops.

Safety

: NiTi resists spontaneous breakage.
Increased angular rotation reduces fracture.
NiTi decreases zipping and stripping.



NITI REAMERS



NICKEL-TITANIUM REAMERS(.02 TAPER)

size	21MM	25MM
#8	NT501-101	NT501-201
#10	NT501-102	NT501-202
#15	NT501-103	NT501-203
#20	NT501-104	NT501-204
#25	NT501-105	NT501-205
#30	NT501-106	NT501-206
#35	NT501-107	NT501-207
#40	NT501-108	NT501-208
#45	NT501-109	NT501-209
#50	NT501-110	NT501-210
#55	NT501-111	NT501-211
#60	NT501-112	NT501-212
#70	NT501-113	NT501-213
#80	NT501-114	NT501-214
#15/40	NT501-191	NT501-291
#45/80	NT501-192	NT501-292

(Pkg. of 6)



NITI K-FILES



NICKEL-TITANIUM K-FILES(.02 TAPER)

size	21MM	25MM
#8	NT502-101	NT502-201
#10	NT502-102	NT502-202
#15	NT502-103	NT502-203
#20	NT502-104	NT502-204
#25	NT502-105	NT502-205
#30	NT502-106	NT502-206
#35	NT502-107	NT502-207
#40	NT502-108	NT502-208
#45	NT502-109	NT502-209
#50	NT502-110	NT502-210
#55	NT502-111	NT502-211
#60	NT502-112	NT502-212
#70	NT502-113	NT502-213
#80	NT502-114	NT502-214
#15/40	NT502-191	NT502-291
#45/80	NT502-192	NT502-292

(Pkg. of 6)



NITI H-FILES



NICKEL-TITANIUM H-FILES(.02 TAPER)

size	21MM	25MM
#8	NT503-101	NT503-201
#10	NT503-102	NT503-202
#15	NT503-103	NT503-203
#20	NT503-104	NT503-204
#25	NT503-105	NT503-205
#30	NT503-106	NT503-206
#35	NT503-107	NT503-207
#40	NT503-108	NT503-208
#45	NT503-109	NT503-209
#50	NT503-110	NT503-210
#55	NT503-111	NT503-211
#60	NT503-112	NT503-212
#70	NT503-113	NT503-213
#80	NT503-114	NT503-214
#15/40	NT503-191	NT503-291
#45/80	NT503-192	NT503-292

(Pkg. of 6)

ENDODONTIC NICKEL TITANIUM HAND FILES/BLISTER PACKAGE



- Blister package for single use
- Increased Flexibility
- Strength
- Convenience
- Efficiency
- Safety

NITI REAMERS
NICKEL-TITANIUM REAMERS(.02 TAPER)

size	21MM	25MM
#8	NTB501-101	NTB501-201
#10	NTB501-102	NTB501-202
#15	NTB501-103	NTB501-203
#20	NTB501-104	NTB501-204
#25	NTB501-105	NTB501-205
#30	NTB501-106	NTB501-206
#35	NTB501-107	NTB501-207
#40	NTB501-108	NTB501-208
#45	NTB501-109	NTB501-209
#50	NTB501-110	NTB501-210
#55	NTB501-111	NTB501-211
#60	NTB501-112	NTB501-212
#70	NTB501-113	NTB501-213
#80	NTB501-114	NTB501-214
#15/40	NTB501-191	NTB501-291
#45/80	NTB501-192	NTB501-292

(Pkg. of 6)

NITI K-FILES
NICKEL-TITANIUM K-FILES(.02 TAPER)

size	21MM	25MM
#8	NTB502-101	NTB502-201
#10	NTB502-102	NTB502-202
#15	NTB502-103	NTB502-203
#20	NTB502-104	NTB502-204
#25	NTB502-105	NTB502-205
#30	NTB502-106	NTB502-206
#35	NTB502-107	NTB502-207
#40	NTB502-108	NTB502-208
#45	NTB502-109	NTB502-209
#50	NTB502-110	NTB502-210
#55	NTB502-111	NTB502-211
#60	NTB502-112	NTB502-212
#70	NTB502-113	NTB502-213
#80	NTB502-114	NTB502-214
#15/40	NTB502-191	NTB502-291
#45/80	NTB502-192	NTB502-292

(Pkg. of 6)

NITI H-FILES
NICKEL-TITANIUM H-FILES(.02 TAPER)

size	21MM	25MM
#8	NTB503-101	NTB503-201
#10	NTB503-102	NTB503-202
#15	NTB503-103	NTB503-203
#20	NTB503-104	NTB503-204
#25	NTB503-105	NTB503-205
#30	NTB503-106	NTB503-206
#35	NTB503-107	NTB503-207
#40	NTB503-108	NTB503-208
#45	NTB503-109	NTB503-209
#50	NTB503-110	NTB503-210
#55	NTB503-111	NTB503-211
#60	NTB503-112	NTB503-212
#70	NTB503-113	NTB503-213
#80	NTB503-114	NTB503-214
#15/40	NTB503-191	NTB503-291
#45/80	NTB503-192	NTB503-292

(Pkg. of 6)

ENDO STOP DISPENSER

• ITEM #701-000

- Contains 120 Assorted Round stops.
- Easy to use.
- Features ruler for measuring.



ENDO STOPS

• ITEM#701-291

- Directional type
- 120 Stops and Assorted colors packed in a plastic box
- These endo stops can not be used with ENDO STOP DISPENSER



ENDODONTIC STAINLESS STEEL HAND FILES

- Root canal cleaning and shaping instruments
- Offers superior cutting efficiency on insertion and withdrawal and procession handling
- Ergonomic color coded handle provides comfort and excellent tactile feedback
- Available in sizes 8 through 80 and 21mm, 25mm & 31mm lengths



REAMERS

- Made from Stainless Steel wire that is ground to triangular shape and twisted

size	21MM	25MM	31MM
# 8	SS501-101	SS501-201	SS501-301
# 10	SS501-102	SS501-202	SS501-302
# 15	SS501-103	SS501-203	SS501-303
# 20	SS501-104	SS501-204	SS501-304
# 25	SS501-105	SS501-205	SS501-305
# 30	SS501-106	SS501-206	SS501-306
# 35	SS501-107	SS501-207	SS501-307
# 40	SS501-108	SS501-208	SS501-308
# 45	SS501-109	SS501-209	SS501-309
# 50	SS501-110	SS501-210	SS501-310
# 55	SS501-111	SS501-211	SS501-311
# 60	SS501-112	SS501-212	SS501-312
# 70	SS501-113	SS501-213	SS501-313
# 80	SS501-114	SS501-214	SS501-314
# 15/40	SS501-191	SS501-291	SS501-391
# 45/80	SS501-192	SS501-292	SS501-392

(Pkg. of 6)

K-FILES

- Designed to clean and shape root canals
- Produced to give operators a smooth tactile sense inside the canal during instrument

size	21MM	25MM	31MM
# 8	SS502-101	SS502-201	SS502-301
# 10	SS502-102	SS502-202	SS502-302
# 15	SS502-103	SS502-203	SS502-303
# 20	SS502-104	SS502-204	SS502-304
# 25	SS502-105	SS502-205	SS502-305
# 30	SS502-106	SS502-206	SS502-306
# 35	SS502-107	SS502-207	SS502-307
# 40	SS502-108	SS502-208	SS502-308
# 45	SS502-109	SS502-209	SS502-309
# 50	SS502-110	SS502-210	SS502-310
# 55	SS502-111	SS502-211	SS502-311
# 60	SS502-112	SS502-212	SS502-312
# 70	SS502-113	SS502-213	SS502-313
# 80	SS502-114	SS502-214	SS502-314
# 15/40	SS502-191	SS502-291	SS502-391
# 45/80	SS502-192	SS502-292	SS502-392

(Pkg. of 6)

H-FILES

- Ground from Stainless Steel wire
- The flute depth increases with the taper, creating a shaper, more flexible instrument


size	21MM	25MM	31MM
# 8	SS503-101	SS503-201	SS503-301
# 10	SS503-102	SS503-202	SS503-302
# 15	SS503-103	SS503-203	SS503-303
# 20	SS503-104	SS503-204	SS503-304
# 25	SS503-105	SS503-205	SS503-305
# 30	SS503-106	SS503-206	SS503-306
# 35	SS503-107	SS503-207	SS503-307
# 40	SS503-108	SS503-208	SS503-308
# 45	SS503-109	SS503-209	SS503-309
# 50	SS503-110	SS503-210	SS503-310
# 55	SS503-111	SS503-211	SS503-311
# 60	SS503-112	SS503-212	SS503-312
# 70	SS503-113	SS503-213	SS503-313
# 80	SS503-114	SS503-214	SS503-314
# 15/40	SS503-191	SS503-291	SS503-391
# 45/80	SS503-192	SS503-292	SS503-392

(Pkg. of 6)

BARBED BROACHES



- Designed to remove pulp tissue
- Put it straight into the root canal, rotation at 180° and pull out



size	21MM	25MM
#20(XXXXF)	504-001	504-201
#25(XXXF)	504-002	504-202
#30(XXF)	504-003	504-203
#35(XF)	504-004	504-204
#40(F)	504-005	504-205
#25/40(XXXF/F)	504-092	504-292

(Pkg. of 6)

SPREADER



- Used for the lateral condensation of root canal fillings.
- Made of Stainless Steel.
- High Resistance to wear and Excellent Cutting Efficiency.
- 6pcs in a blister package

size	21MM	25MM
#15	520-103	520-203
#20	520-104	520-204
#25	520-105	520-205
#30	520-106	520-206
#35	520-107	520-207
#40	520-108	520-208
#15/40	520-191	520-291

PLUGGER



- Used for the vertical condensation of root canal fillings.
- Made of Stainless Steel.
- High Resistance to wear and Excellent Cutting Efficiency.
- 6pcs in a blister package

size	21MM	25MM
#15	521-103	521-203
#20	521-104	521-204
#25	521-105	521-205
#30	521-106	521-206
#35	521-107	521-207
#40	521-108	521-208
#15/40	521-191	521-291

DIA-KONDENSOR


- Having both a nickel-titanium end and a stainless steel end, it can be used in any common endodontic obturation technique including warm and cold condensation.
- Excellent durability.
- Can be supplied 1 piece or a set included each 4 sizes.
- Convenient use with color-coding on the condensers.
- Lightweight and comfort design

size	ITEM No.
#35/70(Green)	551-101
#40/80(Black)	551-102
#50/100(Yellow)	551-103
#60/120(Blue)	551-104
Set (4pcs)	551-191




DIA-ENDORING

- **ITEM #850-001**
134 A
- A convenient ring for Hand & Rotary Files.
- Gauges from 1 to 30mm.
- A simple & practical design for measuring files.
- Gold plated ring prevents color change.
- Dimension: 45 x 25 x 15m/m (W x H x D)



STERI ENDO GUARD W/RULER(ENGINE FILE)

- **ITEM #716-203 Yellow**
- **ITEM #716-206 Blue**
134 P
- Holds 16 NiTi Rotary Files.
- Right or left hand use.
- Available to insert sponge to clean or hold files.
- Diameter: 75m/m, Height: 35m/m




ENDO FINGERING

- **ITEM #790-101**
134 P
- Provides complete control of entire endodontic files.
- Disposable sponges hold files in a convenient position.
- Autoclavable except cups
- Contents :
- **ITEM #790-111 :**
Disposable Cups(48pcs)
- **ITEM #790-121 :**
Disposable Sponges(48pcs)
- Finger Rings(3pcs)
- Cup Rings(2pcs)



FILE MEASURE SET(PLASTIC)

- **ITEM #780-001**
134 P
- Holds 32 Measuring Grades.
- Gauge from 12 to 27.5m/m, Including half measurements.
- Dimensions: 87 x 31 x 20m/m (W x H x D)



DIAFLEX™

IRRIGATION SYRINGE FOR ROOT CANAL



- 134 P
- Flexible soft polypropylene tip syringe for root canal irrigation.
- Autoclavable and chemiclavable for multiple use up to 50 times.
- 1.3ml volume imprinted at intervals of 1mm.
- Provides maximum clinical effect in a minimum of time.
- DiaFlex is compatible with Luer Lock Syringes.


DiaFlex - III

- **ITEM #1003-301**
- Package of 3 color coded syringes, 20 changeable tips and silicone oil

DiaFlex - V


- **ITEM #1003-501**
- Package of 5 color coded syringes, 20 changeable tips and silicone oil

- **ITEM #1001-315**
- Blue 20 disposable tips in a blister package



STERI. IRRIGATION TIPS-SIDE OPEN TYPE

- Disposable Endodontic & Periodontal Irrigation Tips with Safety.
- Unique design with a side-port for superior effectiveness and patient comfort in irrigating root canals and periodontal pockets.
- Pre-sterilized by E.O Gas.
- Use with Lure Lock Syringes.
- Side-Port and Closed-end Irrigation tips.



- PERIODONTAL POCKETS**
- Side-port design make an effective irrigate pockets and closed, rounded tip can reduce the risk of tissue damage.
- ROOT CANALS**
- Side-port dispersal can create an upward flushing motion for the complete root canal preparations and cleanse the gingival sulci fully without damaging tissue.

size	Color	ITEM No.
31G (0.27 x 25mm)	Clear White	951-101
30G (0.30 x 25mm)	Yellow	951-102
28G (0.36 x 25mm)	Blue Green	951-103

(Pkg. of 50)



IRRIGATION TIPS-END OPEN TYPE



< 22G >

< 20G >

- Disposable Endodontic & Periodontal Irrigation tips
- Use with Lure Lock Syringes.
- Thoroughly irrigate root canal preparations.
- Round end tip can prevent the secondary injury during the use.
- Widest selection of tips: 6 sizes from 18G to 27G.

size	Color	ITEM No.
18G (1.2 x 25mm)	Pink	951-201
19G (1.1 x 25mm)	Cream	951-202
20G (0.9 x 25mm)	Yellow	951-203
22G (0.7 x 25mm)	Black	951-204
25G (0.5 x 25mm)	Orange	951-205
27G (0.4 x 25mm)	Medium Grey	951-206

(Pkg. of 100)



MISCELLANEOUS
ITEMS

- Aluminum Endo Storage Boxes
- Protective Shields & Glasses
- Sterilizing Containers & Trays
- Propy Brushes
- Propy Cups
- Bur Blocks & Stands
- Instrument Storage Trays
- Nickel Titanium Arch Wires
- Glass Items
- Dental Wax

SMALL ROUND TYPE A

- ITEM # 803-003 Red
- ITEM # 803-004 Blue



- Has 6 Silicone Color Coded Silicone Cups.
- Ideal for Organizing Gutta Percha Points & Paper Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Diameter: 40m/m, Height: 60m/m



SMALL ROUND TYPE B

- ITEM # 802-003 Red
- ITEM # 802-004 Blue



- Ideal for Organizing Hand Files.
- Holds 20 instruments.
- Removable Insert for Cleaning.
- Diameter: 40m/m, Height: 60m/m



SMALL SQUARE TYPE A

- ITEM #807-003 Red
- ITEM #807-004 Blue



- Has 6 Silicone Color Coded Silicone Cups.
- Ideal for Organizing Gutta Percha Points & Paper Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 50 x 50 x 45m/m (W x H x D)



SMALL SQUARE TYPE B

- ITEM #806-003 Red
- ITEM #806-004 Blue



- Ideal for Organizing Hand Files.
- Holds 30 Endo Instruments.
- Removable Insert for Cleaning.
- Dimensions: 50 x 50 x 45m/m (W x H x D)



MEDIUM OBLONG TYPE A

- ITEM #801-003 Red
- ITEM #801-004 Blue



- A Versatile & Practical Organizer with a Clean Grip for Hand & Rotary File as well as Cups for Gutta Percha & Paper Points.
- Holds 36 Instruments, 3 Silicone Cups, 1 Clean Grip.
- Removable Insert for Cleaning.
- Dimensions: 130 x 50 x 50m/m (W x H x D)



MEDIUM OBLONG TYPE B

- ITEM #809-003 Red
- ITEM #809-004 Blue



- A Versatile & Practical Organizer for Hand & Rotary File as well as Cups for Gutta Percha & Paper Points.
- Holds 60 Instruments, 4 Silicone Cups.
- Removable Insert for Cleaning.
- Dimensions: 130 x 50 x 50m/m (W x H x D)



MEDIUM ROUND

- ITEM # 821-001 Silver



- Ideal for Organizing Hand Files & Gutta Percha Points & Paper Points.
- Holds 64 Instruments, 4 Color Coded Silicone Cups.
- Removable Insert for Cleaning.
- Diameter: 92m/m, Height: 69m/m



MEDIUM SQUARE

- ITEM #804-003 Red
- ITEM #804-004 Blue



- A Versatile & Practical Organizer with a Clean Grip for Hand & Rotary File as well as Cups for Gutta Percha & Paper Points.
- Features a Endo Meter for measuring Files & Points.
- Holds 5 Burs, 30 Instruments, 6 Silicone Cups, 1 Clean Grip, Endo Meter and 1 location for small instruments.
- Removable Insert for Cleaning.
- Dimensions: 120 x 50 x 80m/m (W x H x D)



DIA-DELUXE

- ITEM #811-003 Red



- DiaDent's Ultimate Box for Organizing your Files & Points. Holds 4 Spreaders/Pluggers, 56(21mm) & 48(25mm) Instruments, 6 Gates Glidens, 10 Barbed Broaches, 28 Cups, Section for Cotton.
- Removable Insert for Cleaning.
- Dimensions: 215 x 50 x 100m/m (W x H x D)



LARGE OBLONG TYPE A

- ITEM #812-004 Blue



- A Simple & Effective Organizer with Sections for your Hand Files.
- Holds 180 Endo Instruments, 8 Color Coded Silicone Cups for Points, perforated bottom case.
- Removable Insert for Cleaning.
- Dimensions: 218 x 53 x 105m/m (W x H x D)



SMALL OBLONG-6 COMPT.

- ITEM #805-003 Red
- ITEM #805-004 Blue



- Ideal for Holding Burs & Small Instruments.
- Has 6 Compartments.
- Removable Insert for Cleaning.
- Dimensions: 130 x 15 x 50m/m (W x H x D)



LARGE OBLONG TYPE B

- ITEM #812-005



- An Organizer for Hand & Rotary Files as well as Points.
- Store your Instruments in One Place.
- Holds 357 Endo Instruments, 21 Silicone Cups for Points.
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m (W x H x D)



LARGE OBLONG TYPE C

- ITEM #812-006



- A Very Practical Box. Easy to use inserts individually as needed.
- Holds 206 Endo Instruments, 12 Silicone Cups for Points. Interchangeable Inserts.
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m (W x H x D)



LARGE OBLONG TYPE D

- ITEM #812-007



- Organizer for Hand & Rotary Files.
- Holds 192 Hand Instruments, 80 Rotary Instruments.
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m (W x H x D)



LARGE COMBI

- ITEM # 818-001



- Versatile & Practical. Allows separation of specific instruments.
- Holds 5 individual Inserts(A+B+C+D+E).
- Insert A : 6 cups + 6 cups
- Insert B : 2 sets of 6 cups & Holds 12 Endo Instruments
- Insert C : 4 cups, Holds 48 Endo Instruments
- Insert D : 2 cups, Holds 48 Endo Instruments
- Insert E : Holds 72 Endo Instruments
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m(W x H x D)



MEDIUM COMBI

• ITEM #817-001



- A Compact & Practical File & Point Organizer.
- Holds 3 Individual Inserts(C+D+E).
- Insert C : 4 cups, Holds 48 Endo Instruments
- Insert D : 2 cups, Holds 48 Endo Instruments
- Insert E : Holds 72 Endo Instruments
- Removable Insert for Cleaning.
- Dimensions: 122 x 50 x 100m/m (W x H x D)



SMALL COMBI TYPE B

• ITEM #816-020



- A Compact Organizer for Hand & Rotary Files.
- Hold 24 Endo Instruments, 12 Color Coded Silicone Cups for Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE D

• ITEM #816-040



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments, 2 Color Coded Silicone Cups for Points.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE F

• ITEM #816-060



- A Compact Organizer for Gutta Percha & Paper Points.
- Holds 12 Color Coded Silicone Cups.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE A

• ITEM # 816-010



- A Compact organizer for Gutta Percha & Paper Points.
- Holds 12 Color Coded Silicone Cups.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE C

• ITEM #816-030



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments, 4 Color Coded Silicone Cups for Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE E

• ITEM #816-050



- A Compact Organizer for Hand & Rotary Files.
- Holds 72 Endo Instruments.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE G

• ITEM #816-070



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments Medium Sized Perforations.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



SMALL COMBI TYPE H

• ITEM #816-080



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments Large Sized Perforations.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



ENDO GRIP(MINI)

• ITEM #833-001



- Commonly Used for Cleaning Files.
- Holds sponge or gauze.
- Available in Blue.
- Diameter: 40m/m, Height: 35m/m



ENDO METER SQUARE

• ITEM # 835-001



- Holds & Measures 32 Files or Other Endo Instruments.
- Gauge from 15 to 30.5m/m, including half measurement.
- 2 Cups for Files, 1 location for Endo Stops.
- Dimensions: 125 x 40 x 40m/m (W x H x D)



ENDO METER

• ITEM #835-002 Type A- Silver
• ITEM #835-003 Type B- Red



- Holds & Measures 10 Files.
- Gauges from 1 to 35mm.
- Features a locking spring mechanism inside for easy use.
- Diameter: 35m/m, Height: 48m/m



Type A



Type B

SMALL COMBI TYPE I

• ITEM #816-090



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments Medium & Large Sized Perforations.
- Removable Insert for Cleaning.
- Dimensions: 105 x 52 x 44m/m (W x H x D)



ENDO CLEANER

• ITEM #834-001 Open
• ITEM #834-002 Closed



- Commonly Used for Cleaning and sterilizing Files.
- Top unscrews for easy replacement of Gauze.
- Diameter: 40m/m, Height: 30m/m



LARGE OBLONG EMPTY BOX

• ITEM #820-003 Blue



- An all purpose storage container.
- Ideal for holding accessories.
- Dimensions: 210 x 25 x 50m/m (W x H x D)



DIA-100FS/DIA-100SS

- Optical grade protective shields and glasses for the face and eyes from flying debris or substances.
- Lightweight and very comfortable.
- Shields can be quickly positioned and fits over eyeglasses
- Easy to clean, washable

Dia-100FS

- Full face protection shield
- 1 Frame & 5 shields

1005-5010	DIA-100FS	REFILL PACKAGE	5 shields only
1005-5011	DIA-100FS	REGULAR PACKAGE	CRYSTAL
1005-5012	DIA-100FS	REGULAR PACKAGE	WHITE
1005-5014	DIA-100FS	REGULAR PACKAGE	BLUE

DIA-100SS

- Eye protection glasses
- 1 Frame & 5 shields

1005-5020	DIA-100SS	REFILL PACKAGE	5 shields only
1005-5021	DIA-100SS	REGULAR PACKAGE	CRYSTAL
1005-5022	DIA-100SS	REGULAR PACKAGE	WHITE
1005-5024	DIA-100SS	REGULAR PACKAGE	BLUE



DIA-100FS




DIA-100SS

DIA-300I

- **ITEM #1005-1041**
- Protection Glasses. Protect your eyes from flying debris during clinical or laboratory procedures.
- High optical quality lenses with chemical resistance properties.
- Comfort fit at the nose bridge at the temple.
- Lens coating(PC lenses, Anti-fog, Anti-UV, Anti-static)

Product Features



Vertical adjustments

Horizontal adjustments

Ultra light weight



Anti-fog
Anti-UV
Anti-static



DIA-200B

- **ITEM #1005-1031**
- Multi Purpose Goggles. Used among dentists and surgeons for eyes protection from flying debris. Preventing eye injuries and contamination from fine dust.
- Spherical lens design for superior optics.
- Comfortable lightweight fit.
- Fits over Eyeglasses.
- Designed for coolness and ventilation.
- Anti-fog coating, Anti-static.

Product Features



Superior cushioning and ventilation

Good air circulation

Easily used with correction glasses



Anti-fog
Anti-static



GERMICIDE TRAY

- **ITEM #711-001**



- Used for Cold Sterilization.
- Easy & Safe to handle.
- One Large & One Small compartment.
- Easy loading and unloading.
- Dimensions: 263 x 72 x 115 m/m (W x H x D)



CONTAINER(STERILIZING)



- **ITEM #824-001 Type A**
- Holds 10 Endo Instruments.
- Dimensions : 185 x 35 x 145 m/m (W x H x D)



Type A

- **ITEM #824-002 Type B**
- Holds 16 Endo Instruments.
- Dimensions : 185 x 35 x 286 m/m (W x H x D)



Type B

- Well Built, Rigid, Sterilizing Container.
- Features Locking System, Fits many sterilizing units.

STERILIZING DRUM

- **ITEM # 828-002 Large**



- Diameter: 115m/m, Height: 100m/m

Autoclaving burs & smaller struments
Drainage



Large

- **ITEM #814-001 Small**



- Diameter: 56m/m, Height: 47m/m
- Color: Silver

Separates into 3 pieces for easy Cleaning.
Equipped with locking system.



Small

FLAT TRAY(STERILIZING) TYPE A

- **ITEM #827-001 Blue**



- Sterilizing Tray.
- Ideal for Operator Instruments.
- Anodized Aluminum.
- Dimensions : 200 x 20 x 105m/m (W x H x D)



FLAT TRAY(STERILIZING) TYPE B

- **ITEM #827-002 Red**



- Sterilizing Tray.
- Ideal for Operator Instruments.
- Anodized Aluminum.
- Dimensions : 180 x 10 x 100m/m (W x H x D)



GERMICIDE TRAY TYPE B

- **ITEM #711-106 Blue**



- Efficient high-volume germicide tray designed for maximum use of chemicals to clean/disinfect/sterilize instruments
- Dimensions : 350 x 150 x 250 m/m (W x H x D)



CONTAINER(STORAGE)



- **ITEM #829-001 Type A**
- Dimensions: 255 x 50 x 216m/m (W x H x D)



Type A

- **ITEM #829-002 Type B**
- Dimensions: 255 x 50 x 108m/m (W x H x D)




Type B

- Multi Purpose Box.
- Use as an organizer or for storage.
- Bur Holders fit well in this container.

WASTE CONTAINER

SCALING BRUSH LATCH TYPE

- **ITEM #941-001**
- Packed 72 brushes.
- Regular cup style.
- For polishing of tooth surface after scaling.
- Useful for wide space.
- Diameter: 4.49mm



SCALING POLISHING BRUSH LATCH TYPE

- **ITEM #941-103**
- Packed 36 brushes.
- Small cup style.
- Useful for narrow space.
- Diameter : 2.27mm




CANDLE BRUSH

- **ITEM #941-104 :**
Latch Type, Packed 24 brushes
- **ITEM #941-105 :**
Screw Type, Packed 24 brushes and 10 mandrels
- For removing stains from occlusal pits and fissures and for preparing tooth surface prior to sealant and orthodontic application.
- Diameter : Latch Type : 4.50mm / Screw Type : 5.00mm



POST BRUSH

- **ITEM #941-101 :**
Latch Type, Packed 30 brushes
- **ITEM #941-102 :**
Screw Type, Packed 30 brushes
- For quick cleaning root canal preparation, under bridge.
- Available to insert sterilizing liquid during turn it left in Contra Angle.
- Available to flow out the stains like cement, liquid during turn it right.
- Recommended rotation: 300-800rpm.
- Single use.
- Diameter :
Latch Type : 1.20mm / Screw Type : 1.03mm



PROPHY BRUSH

- **ITEM #950-101 :** Latch Type
- **ITEM #950-102 :** Screw Type
- Packed 144 brushes
- Diameter
Latch Type : 4.48mm
Screw Type : 4.48mm



Latch Type (PB - 001)



Screw Type (PB - 002)

PROPHY BRUSH POINTED HEAD

- **ITEM #950-111 :** Latch Type
- **ITEM #950-112 :** Screw Type
- Packed 144 brushes
- Diameter
Latch Type : 4.48mm
Screw Type : 4.48mm



Latch Type (PB - 003)



Screw Type (PB - 004)

PROPHY BRUSH LARGE HEAD LATCH TYPE

- **ITEM #950-121**
- Packed 144 brushes
- Diameter : 5.98mm



(PB - 005)

PROPHY BRUSH LARGE POINTED HEAD LATCH TYPE

- **ITEM #950-131**
- Packed 144 brushes
- Diameter : 5.98mm



(PB - 006)

SILICONE CUP FOR POLISHING LATCH TYPE

- **ITEM #941-002**
- Soft, white cups
- For polishing the tooth surface, achieves a high shine quickly.
- Use conveniently after insert in Contra Angle
- Pkg. of 72



PROPHY RUBBER CUP RIBBED & WEBBED

- **ITEM #950-201 :** Latch Type
- **ITEM #950-202 :** Screw Type
- Packed 144 cups



Latch Type (RC - 001)



Screw Type (RC - 002)

PROPHY RUBBER CUP WEBBED

- **ITEM #950-211 :** Latch Type
- **ITEM #950-212 :** Screw Type
- Packed 144 cups



Latch Type (RC - 003)



Screw Type (RC - 004)

PROPHY RUBBER CUP RIBBED

- **ITEM #950-221 :** Latch Type
- **ITEM #950-222 :** Screw Type
- Packed 144 cups



Latch Type (RC - 005)



Screw Type (RC - 006)

PROPHY RUBBER CUP UMBRELLA SHAPE

- **ITEM #950-231 :** Latch Type
- **ITEM #950-232 :** Screw Type
- Packed 144 cups



Latch Type (RC - 007)



Screw Type (RC - 008)

PROPHY RUBBER CUP SNAP-ON

- **ITEM #950-251 :** Ribbed
- **ITEM #950-252 :** Webbed
- Packed 144 cups



Ribbed (RC - 010)



Webbed (RC - 011)

PROPHY RUBBER CUP WEBBED WITH WASHER LATCH TYPE

- **ITEM #950-241**
- Packed 144 cups



(RC - 009)

PROPHY RUBBER CUP SNAP-ON RIBBED & WEBBED

- **ITEM #950-253**
- Packed 144 cups



(RC - 012)

STERI BUR BLOCK W/COVER

Type A

• ITEM #837-001

- Holds 24 Friction Grip(high speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 × 45 × 44m/m (W × H × D)



Type B

• ITEM #837-002

- Holds 24 Right Angle(low speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 × 45 × 44m/m(W × H × D)



Type C

• ITEM #837-003

- Holds 12 Friction Grip(high speed) Burs & 12 Right Angle(low speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 × 45 × 44m/m(W × H × D)



Type D

• ITEM #837-004

- Holds 36 Friction Grip(high speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 × 45 × 50m/m (W × H × D)



Type E

• ITEM #837-005

- Holds 36 Right Angle(low speed) Burs.
- Both plastic cover and bur block are Autoclavable.
- Dimensions : 60 × 45 × 50m/m (W × H × D)




MINI BUR STAND

Type A



- Silver ITEM #838-101
- Red ITEM #838-103
- Blue ITEM #838-104

- Holds 12 Friction Grip(high speed) Burs
- Durable Aluminum
- Dimensions: 45 × 32 × 25m/m (W × H × D)



Type B



- Silver ITEM #838-201
- Red ITEM #838-203
- Blue ITEM #838-204

- Holds 12 Right Angle(low speed) Burs
- Durable Aluminum
- Dimensions: 45 × 32 × 25m/m (W × H × D)



Type C



- Silver ITEM #838-301
- Gold ITEM #838-302
- Red ITEM #838-303
- Blue ITEM #838-304

- Holds 9 Friction Grip(high speed) Burs & 3 Right Angle(low speed) Burs
- Durable Aluminum
- Dimensions: 45 × 32 × 25 m/m (W × H × D)



※ All of Mini Bur Stands can be supplied with the private brand after laser markings on each stand.

GERMICIDE ROUND BUR BOX

• ITEM #915-001



- Contains perforated insert for drainage.
- Practical & Durable.
- 1 Basin bowl/ 1 Perforated insert/1 Cover.
- Diameter: 75m/m, Height: 35m/m



DIVIDED TRAY(SIZE B)

- ITEM #715-403 Yellow
- ITEM #715-406 Blue



- Sturdy, non-porous/stain-resistant polypropylene plastic trays
- Available in divided style
- Dimensions : 338 × 23 × 243m/m (W × H × D)



INSTRUMENT TRAY TYPE A (Perforated)

• ITEM #826-001



- Holds 10 Instruments
- Ideal for Instrument Sterilization.
- Perforations Facilitate Sterilization & Drainage.
- Dimensions: 190 × 35 × 150m/m (W × H × D)




INSTRUMENT TRAY TYPE B (Perforated)

• ITEM #826-002



- Holds 16 Instruments
- Ideal for Instrument Sterilization.
- Perforations Facilitate Sterilization & Drainage.
- Dimensions: 280 × 35 × 180m/m (W × H × D)



COMPARTMENT BOX TYPE A

• ITEM #822-001



- Holds 16 basins.
- Small Item Organizer.
- Dimensions: 230 × 25 × 110m/m (W × H × D)



STAINLESS STEEL TRAY TYPE A

• ITEM #914-001



- A larger Tray for all purposes and storing.
- Dimensions: 180 × 40 × 135m/m (W × H × D)



INSTRUMENT TRAY TYPE A(Non-Perforated)

• ITEM #825-001



- Holds 10 Instruments
- Ideal for Instrument Storage.
- Dimensions: 190 × 35 × 150m/m (W × H × D)



INSTRUMENT TRAY TYPE B (Non-Perforated)

• ITEM #825-002



- Holds 16 Instruments
- Ideal for Instrument Sterilization.
- Dimensions: 280 × 35 × 180m/m (W × H × D)



COMPARTMENT BOX TYPE B

• ITEM #822-002




- Holds 12 basins.
- Compartments vary in Size.
- Dimensions: 230 × 25 × 110m/m (W × H × D)



NICKEL TITANIUM ARCH WIRES

- Last up to 5 times more than Stainless steel
- Larger working range than Stainless steel
- Less distortion and more shape memory than Stainless steel
- Easy to arrange in the bracket slots



SUPER ELASTIC ARCH WIRES

- Made of the highest quality Nickel-Titanium wire
- Highly polished ultra smooth finishing
- Allow a high degree of deflection
- Comes in Round(for class 1, mild cases) and Rectangular (for more stiffness and three dimensional control) Wires.

Round Type: Pkg. of 10

NT311-121	.012 Upper	SE-D-12-U
NT311-122	.012 Lower	SE-D-12-L
NT311-141	.014 Upper	SE-D-14-U
NT311-142	.014 Lower	SE-D-14-L
NT311-161	.016 Upper	SE-D-16-U
NT311-162	.016 Lower	SE-D-16-L
NT311-181	.018 Upper	SE-D-18-U
NT311-182	.018 Lower	SE-D-18-L
NT311-201	.020 Upper	SE-D-20-U
NT311-202	.020 Lower	SE-D-20-L

Rectangular Type: Pkg. of 10

NT321-661	.016 x .016 Upper	SE-R-1616-U
NT321-662	.016 x .016 Lower	SE-R-1616-L
NT321-621	.016 x .022 Upper	SE-R-1622-U
NT321-622	.016 x .022 Lower	SE-R-1622-L
NT321-721	.017 x .022 Upper	SE-R-1722-U
NT321-722	.017 x .022 Lower	SE-R-1722-L
NT321-751	.017 x .025 Upper	SE-R-1725-U
NT321-752	.017 x .025 Lower	SE-R-1725-L
NT321-821	.018 x .022 Upper	SE-R-1822-U
NT321-822	.018 x .022 Lower	SE-R-1822-L
NT321-851	.018 x .025 Upper	SE-R-1825-U
NT321-852	.018 x .025 Lower	SE-R-1825-L
NT321-951	.019 x .025 Upper	SE-R-1925-U
NT321-952	.019 x .025 Lower	SE-R-1925-L
NT321-151	.021 x .025 Upper	SE-R-2125-U
NT321-152	.021 x .025 Lower	SE-R-2125-L



REVERSE CURVE

Round Type: Pkg. of 10

NT313-121	.012 Upper	RC-D-12-U
NT313-122	.012 Lower	RC-D-12-L
NT313-141	.014 Upper	RC-D-14-U
NT313-142	.014 Lower	RC-D-14-L
NT313-161	.016 Upper	RC-D-16-U
NT313-162	.016 Lower	RC-D-16-L
NT313-181	.018 Upper	RC-D-18-U
NT313-182	.018 Lower	RC-D-18-L
NT313-201	.020 Upper	RC-D-20-U
NT313-202	.020 Lower	RC-D-20-L

Rectangular Type: Pkg. of 10

NT323-661	.016 x .016 Upper	RC-R-1616-U
NT323-662	.016 x .016 Lower	RC-R-1616-L
NT323-621	.016 x .022 Upper	RC-R-1622-U
NT323-622	.016 x .022 Lower	RC-R-1622-L
NT323-721	.017 x .022 Upper	RC-R-1722-U
NT323-722	.017 x .022 Lower	RC-R-1722-L
NT323-751	.017 x .025 Upper	RC-R-1725-U
NT323-752	.017 x .025 Lower	RC-R-1725-L
NT323-821	.018 x .022 Upper	RC-R-1822-U
NT323-822	.018 x .022 Lower	RC-R-1822-L
NT323-851	.018 x .025 Upper	RC-R-1825-U
NT323-852	.018 x .025 Lower	RC-R-1825-L
NT323-951	.019 x .025 Upper	RC-R-1925-U
NT323-952	.019 x .025 Lower	RC-R-1925-L
NT323-151	.021 x .025 Upper	RC-R-2125-U
NT323-152	.021 x .025 Lower	RC-R-2125-L

SUPER ELASTIC OPEN SPRING


Pkg. of 2

NT331-010	.010	SE-S-10-O
NT331-012	.012	SE-S-12-O

SUPER ELASTIC CLOSE SPRING

Pkg. of 10

NT332-010	.010	SE-S-10-C
NT332-012	.012	SE-S-12-C



THERMAL ACTIVATED ARCH WIRES

- Achieve gentle torque force at the early stage of treatment
- Achieve optimum activation at 37.5°C (Mouth Temperature) and then the force level is increased
- Easy Manipulation

Round Type: Pkg. of 10

NT312-121	.012 Upper	SM-D-12-U
NT312-122	.012 Lower	SM-D-12-L
NT312-141	.014 Upper	SM-D-14-U
NT312-142	.014 Lower	SM-D-14-L
NT312-161	.016 Upper	SM-D-16-U
NT312-162	.016 Lower	SM-D-16-L
NT312-181	.018 Upper	SM-D-18-U
NT312-182	.018 Lower	SM-D-18-L
NT312-201	.020 Upper	SM-D-20-U
NT312-202	.020 Lower	SM-D-20-L

Rectangular Type: Pkg. of 10

NT322-661	.016 x .016 Upper	SM-R-1616-U
NT322-662	.016 x .016 Lower	SM-R-1616-L
NT322-621	.016 x .022 Upper	SM-R-1622-U
NT322-622	.016 x .022 Lower	SM-R-1622-L
NT322-721	.017 x .022 Upper	SM-R-1722-U
NT322-722	.017 x .022 Lower	SM-R-1722-L
NT322-751	.017 x .025 Upper	SM-R-1725-U
NT322-752	.017 x .025 Lower	SM-R-1725-L
NT322-821	.018 x .022 Upper	SM-R-1822-U
NT322-822	.018 x .022 Lower	SM-R-1822-L
NT322-851	.018 x .025 Upper	SM-R-1825-U
NT322-852	.018 x .025 Lower	SM-R-1825-L
NT322-951	.019 x .025 Upper	SM-R-1925-U
NT322-952	.019 x .025 Lower	SM-R-1925-L
NT322-151	.021 x .025 Upper	SM-R-2125-U
NT322-152	.021 x .025 Lower	SM-R-2125-L

GLASS DAPPEN DISHES

- ITEM #601-001 Clear
- ITEM #601-002 Amber
- ITEM #601-003 Blue
- ITEM #601-004 Green



- Standard Dappen Dish.
- Do not Autoclave.
- Easy to clean.
- Sold Individually or Box of 12pcs.



MEDICAMENT BOTTLES



ROUND TYPE

- ITEM #604-001 Clear
- ITEM #604-002 Amber
- ITEM #604-003 Blue
- ITEM #604-004 Green

SQUARE TYPE

- ITEM #605-001 Clear
- ITEM #605-002 Amber
- ITEM #605-003 Blue
- ITEM #605-004 Green

- Locking lids for good sealing capabilities.
- Available in 4 colors for coding of chemicals.
- Easy to use and store.
- High quality glass, economically priced.
- 10ml volume
- Do not Autoclave.



MIXING SLAB



- ITEM #603-005 LT
- Dimensions : 150 x 15 x 80m/m (W x H x D)
- ITEM #603-004 LM
- Dimensions : 150 x 10 x 75m/m (W x H x D)
- ITEM #603-003 M
- Dimensions : 122 x 5 x 75m/m (W x H x D)
- ITEM #603-002 L
- Dimensions : 112 x 6 x 76m/m (W x H x D)
- ITEM #603-001 S
- Dimensions : 76 x 6 x 65m/m (W x H x D)
- High quality, competitively priced Mixing Slabs available in 5 sizes.
- Do not Autoclave.



CHAIR SIDE TABLE

- ITEM #712-001 White



- Comes with attachments.
- Can be adapted to several types of chairs.
- Hole Diameter : available from 25 to 50m/m.
- Easy access to instruments.
- Dimensions : 390 x 30 x 235m/m (W x H x D)



GLASS ENDO ORGANIZER

- ITEM #602-000



- Shallow compartment walls for liquid sterilization.
- 6 inch diameter dish with 12 shallow wall compartments.
- Glass lid with depressions for chemicals.
- Do not Autoclave.



MEDICAMENT BOTTLE STAND

- ITEM #606-001



- Fits 5 Round or Square Medicament Bottles.
- Opaque high quality glass.
- Do not Autoclave.
- Color : Milky White
- Dimensions : 215 x 27 x 55m/m (W x H x D)



COTTON ROLL DISPENSER

- ITEM #770-001 White



- Holds approximately 60 Cotton Rolls.
- Features pull out drawer for easy cotton dispensing.
- Lid flips back for refilling.
- Non-slip bottom.
- Adjustable width for Cotton Roll in 20, 30 and 40m/m.
- Dimensions : 170 x 101 x 52m/m (W x H x D)



TISSUE HOLDER

- ITEM #712-101 White



- Can be adapted to several types of chairs.
- Tissue box holder fits on unit chair pole.
- Hole Diameter : available from 25 to 50m/m.
- Dimensions : 245 x 54 x 125m/m (W x H x D)



STICKY WAX

- **ITEM # 940-010 Beige**
- Made of superior resins and wax.
- Used for the temporary joint in the denture repair works.
- Easy to scrape off the denture after the plaster case is made.
- Diameter: 7m/m, Height: 80m/m
- 100g/box



UTILITY WAX

- **ITEM # 940-020 White**
- **ITEM # 940-021 Red**
- Used to build a bank along the tray.
- Has the superior flexibility
- Diameter: 5m/m, Height: 315m/m
- 250g/box



INLAY BLOCK WAX

- **ITEM # 940-030 Green**
- Used for Inlay casting.
- Low congelation and contraction, and adequate cohesion and solidity levels make it possible to produce a more precise waxing pattern.
- Dimensions: 71 x 11 x 78m/m (W x H x D)
- 450g/box



MODELING WAX

- **ITEM # 940-040 Pink**
- Used for the denture works.
- Provided the sufficient stability for the base plate.
- The solidity remains equal both in the oral and room temperature.
- Dimensions: 170 x 1.4 x 100m/m (W x H x D)
- 500g/box



WAX SEPARATOR

- **ITEM # 940-050 Yellow**
- Used for wax separation.
- Quick drying and not to stick to the denture surface, easily can be separated after some wax-ups.
- 2 bottles of 60ml and 1 bottle of 15ml
- This item should be packed as "Dangerous Goods Packaging" for the shipment.



DIE HARDENER

- **ITEM # 940-070 Blue**
- Makes the surface of abutment stone harden to be consolidated and coated and prevent to be abrasive.
- 2 bottles of 60ml and 1 bottle of 15ml
- This item should be packed as "Dangerous Goods Packaging" for the shipment.



COLOR SPACER

- **ITEM #940-060 SILVER**
- **ITEM #940-061 GOLD**
- **ITEM #940-062 RED**
- One round of application generates space of 12 and 13µm.
- Before use, shake the bottle of Color Spacer.
- When the solution gets stiff, add a few drops of Diluent.
- 1 bottle of 60ml, 1 bottle of 15ml and 1 bottle(Diluent) of 60ml
- This item should be packed as "Dangerous Goods Packaging" for the shipment.



SILVER



GOLD



RED

X-RAY FILM HOLDER

- **ITEM #750-001**
- Used for Anterior and Posterior X-Ray.
- Angled for patient comfort in Anterior x-rays.



RETRACTOR(TYPE A)

- **ITEM #740-001 Adult**
- **ITEM #740-002 Child**

- Soft Flexible Plastic
- Pkg. of 2



ID-RING

- **ITEM #942-101 White**
- **ITEM #942-102 Black**
- **ITEM #942-103 Brown**
- **ITEM #942-104 Green**
- **ITEM #942-105 Orange**
- **ITEM #942-106 Red**
- **ITEM #942-107 Gray**
- **ITEM #942-108 Yellow**
- **ITEM #942-109 Pink**
- **ITEM #942-110 Light Brown**
- **ITEM #942-191 Assorted**
- Silicone rubber
- Standard size(inner diameter 3.6mm)
- Pkg. of 120



CHEEK RETRACTOR(TYPE B)

- **ITEM #741-001 Adult**
- **ITEM #741-002 Child**

- Soft Flexible Plastic
- Expands the mouth widely
- Pkg. of 2



ENDO TRAINING BLOC

- **ITEM #943-100**
- Clear acrylic blocs.
- Ideal for training bloc for Warm Gutta Percha Obturation techniques.
- Pkg. of 10



DOWEL PIN

- **ITEM #907-001: Twin pin**
- **ITEM #907-002: Single pin**

Twin pin with metal sleeve & Rubber cap

- Rubber cap saves working time.
- Ruber cap facilitates work.
- Exact fit into preparations and no distortion.
- Pkg.of 1,000

